



Blader Retractors
Urinary Instruments & Trocars




15

PAGE 204 TO 208


www.humcorp.com

15-104
57x25mm

15-105
82x25mm




15-101
JUDD-MASSON
Bladder-and suprapubic
retractor



15-102
MILLIN
Center blade for (15-101)
retro pubic prostatectomy




15-103
MILLIN
Bladder retractor
15-103A Spreader only
15-104 Blades 57x25mm
15-105 Blades 82x25mm
15-106 Center blade 105x45mm




15-107
LEGUEU
Bladder retractor
Center blade
80 x 31 mm




15-108
LEGUEU
Bladder spatula
260 mm



15-109
KOCHER
Bladder spatula
235 mm




15-110
MATHIEU
Foreign body forceps
Length of shaft
190 mm




15-111
MILLIN
Lobe grasping
forceps
235 mm

Atraumata



15-112
Ureter forceps
240 mm

Atraumata



15-113
MILLIN
Lobe grasping
forceps
245 mm



15-114
MILLIN
Bladder neck spreader
285 mm

15-115
GUYON
Catheter guide
(introducer)
365 mm

1/2

**Charr.
(French)**

15-116	6
15-117	7
15-118	8
15-119	9
15-120	10
15-121	11
15-122	12
15-123	13
15-124	14
15-125	15
15-126	16
15-127	18
15-128	20

1/2

**Charr.
(French)**

15-129	8
15-130	9
15-131	10
15-132	11
15-133	12
15-134	13
15-135	14
15-136	15
15-137	16
15-138	17
15-139	18
15-140	19
15-141	20
15-142	21
15-143	22
15-144	23
15-145	24
15-146	25
15-147	26
15-148	27
15-149	28
15-150	29
15-151	30

1/2

15-152	Ben. 24/Charr. 12
15-153	Ben. 25
15-154	Ben. 26/Charr. 13
15-155	Ben. 27
15-156	Ben. 28/Charr. 14
15-157	Ben. 29
15-158	Ben. 30/Charr. 15
15-159	Ben. 31
15-160	Ben. 32/Charr. 16
15-161	Ben. 33
15-162	Ben. 34/Charr. 17
15-163	Ben. 35
15-164	Ben. 36/Charr. 18
15-165	Ben. 37
15-166	Ben. 38/Charr. 19
15-167	Ben. 39
15-168	Ben. 40/Charr. 20
15-169	Ben. 41
15-170	Ben. 42/Charr. 21
15-171	Ben. 43
15-172	Ben. 44/Charr. 22
15-173	Ben. 45
15-174	Ben. 46/Charr. 23
15-175	Ben. 47
15-176	Ben. 48/Charr. 24
15-177	Ben. 49
15-178	Ben. 50/Charr. 25
15-179	Ben. 51
15-180	Ben. 52/Charr. 26
15-181	Ben. 53
15-182	Ben. 54/Charr. 27
15-183	Ben. 55
15-184	Ben. 56/Charr. 28
15-185	Ben. 57
15-186	Ben. 58/Charr. 29
15-187	Ben. 59
15-188	Ben. 60/Charr. 30

1/2


VEN BUREN
15-115-8 to 15-115-36
Set of 29

15-116 - 128
Female catheters,
metal
160 mm

15-129 - 151
DITTEL
Urethral bougies
and dilating sounds
345 mm

15-152 - 188
GUYON
Urethral bougies
and dilating sounds,
without filiform guides,
with thread in tip
255 mm

Genito-Urinary Instruments & Circumcision Clamps



**GUYON
15-189**


Set of 29; Fig. 8 to Fig. 36

Fr.	Fr.	Fr.
15-189-08	15-189-19	15-189-30
15-189-09	15-189-20	15-189-31
15-189-10	15-189-21	15-189-32
15-189-11	15-189-22	15-189-33
15-189-12	15-189-23	15-189-34
15-189-13	15-189-24	15-189-35
15-189-14	15-189-25	15-189-36
15-189-15	15-189-26	
15-189-16	15-189-27	
15-189-17	15-189-28	
15-189-18	15-189-29	

Assisted by two haemostats or ligatures the foreskin is pushed over the bellshaped protection - rod and the glans is placed into it.
By now, the wing - nut is tightened as close as possible, so the foreskin in annular compressed. The foreskin is cut of with a scalpel, right nest to the compression ring. The circumcision clamp remains a few minutes to avoid bleeding. In case that any large vein, specially the frenular artery, is cut through, a ligature must be set.
The advantage of the method described is the cosmetically favourable scar which is left after rejection of the small necrotized portion.



**MOGEN
15-192A
7 cm**




15-190 - 192

WINKELMANN


15-190	12mm Ø
15-191	16mm Ø
15-192	20mm Ø

Circumcision Clamps



15-190 A	8 mm Ø
15-190 B	11 mm Ø
15-190 C	13 mm Ø
15-190 D	16 mm Ø
15-190 E	21 mm Ø
15-190 F	26 mm Ø
15-190 G	32 mm Ø

Luer Lock



15-193
KLIMA ROSEGGER
35 X 1,8 mm
Sternal puncture needle

15-194
Universal trocar
with 4 tips



15-195 - 199
NELSON

	Inside diam	Outside diam
15-195	2,2 mm	3,0 mm
15-196	3,0 mm	4,0 mm
15-197	4,0 mm	5,0 mm
15-198	5,0 mm	6,0 mm
15-199	5,5 mm	6,5 mm


for hydrocele and abdominal
puncture

15-200 - 203

	Inside diam	Outside diam
15-200	2,2 mm	3,0 mm
15-201	3,0 mm	4,0 mm
15-202	4,0 mm	5,0 mm
15-203	5,0 mm	6,0 mm

for abdominal puncture,
with stop-cock and lateral
drain


15-204 100 mm
15-205 110 mm



15-204 - 205
BUELAU

	Inside diam	Outside diam
15-204	6,5 mm	8,0 mm
15-205	9,0 mm	10,0 mm

Empyema trocar



15-206 - 209
OCHSNER

	Inside diam	Outside diam
15-206	3,0 mm	4,0 mm
15-207	3,6 mm	4,6 mm
15-208	4,3 mm	5,3 mm
15-209	5,0 mm	6,0 mm

Gall bladder trocars

Record Cone 1



15-210
PIERCE
Antrum trocar
ø 4,0 mm

Luer Lock Cone



15-211
COAKLEY
Antrum trocar
30 mm