

Biopsy Instruments




19

PAGE 237 TO 245

www.humcorp.com

Suction Curettes, Uterine Biopsy Curettes



19-103A
Without basket

19-103B
With basket

Malleable



1/4

Malleable



1/4

MOD.MUENCHEN
19-102A
27 cm
4.8 mm ϕ

Stumpf blunt



1/4

Stumpf blunt



1/4



1/4



KEVORKIAN
30 cm / 12"



1/4

Luer-Lock

Luer-Lock


237

19-101
ASCH
Uterine
Secretion scoop
200 mm


19-102
UTERINE
Secretion scoop
270 mm

19-103
NOVAK-SCHOECKAERT
Endometrial biopsy
curette
240 mm

19-104
RANDALL
Endometrial
biopsy curette
240 mm




Biopsy Forceps



19-113 - 115 SEIDL DBP
Biopsy forceps, bayonet shaped for portio, vagina and vulva 240mm

- 2 Different shapes of the jaw are available
- a) Gothic shape 2.0 resp 2.4mm wide (causes less bleeding)
 - b) Rectangular shape, 2.4mm wide



Uterine Scoops



Table of Recamier		Sharp		Blunt
Fig	Width	Rigid	Malleable	Rigid
00	4.5mm	19.130		
0	7.5mm	19.131		19.142
1	8.5mm	19.132	19.153	19.143
2	10.5mm	19.133	19.154	19.144
3	12.0mm	19.134	19.155	19.145
4	14.0mm	19.135	19.156	19.146
5	15.0mm	19.136	19.157	19.147
6	16.5mm	19.137	19.158	19.148
7	18.0mm	19.138	19.159	19.149
8	19.5mm	19.139	19.160	19.150
9	21.5mm	19.140		19.151
10	23.5mm	19.141		19.152
19-130 - 19-141		Sharp - Rigid		
19-153 - 19-160		Sharp - Malleable		
19-142 - 15-152		Blunt - Rigid		





Fig:2 10.5mm
Fig:3 12mm
Fig:4 14mm
Fig:5 15mm





Fig:6 16.5mm
Fig:7 18mm
Fig:8 19.5mm



Fig:9 21.5mm
Fig:10 23.5mm



19-128 - 128E
SIMON
max.240 mm


19-129 - 129A
SCHRODER
max.310 mm

19-130 - 160
RECAMIER
max.300 mm

Fig	Width	Sharp		Blunt
		Rigid	Malleable	Rigid
1	7.5mm	19.161	19.167	19.173
2	8.5mm	19.162	19.168	19.174
3	10.5mm	19.163	19.169	19.175
4	12.0mm	19.164	19.170	19.176
5	14.0mm	19.165	19.171	19.177
6	15.0mm	19.166	19.172	19.178
19-161 - 19-166		Sharp - Rigid		
19-167 - 19-172		Sharp - Malleable		
19-173 - 19-178		Blunt - Rigid		

Fig	Width	Blunt - Malleable
0	16.5mm	19.179
1	18.5mm	19.180
2	19.5mm	19.181
3	21.5mm	19.182
4	23.5mm	19.183
5	26.0mm	19.184
6	28.5mm	19.185

BUMM



Placenta-and Ovum Curettes



Fig	Width	Sharp Malleable	Blunt Malleable
1	23.5mm	19.190	19.192
2	26.0mm		19.193
3	28.5mm	19.191	19.194



19-187
COLLIN
300 mm



19-188 - 189
RECAMIER
blunt-rigid
max. 360 mm




Fig:1
19.188

Fig:2
19.189



19-190 - 194
BUMM
max. 310 mm




Fig:1

Fig:2

Fig:3

sharp-malleable
blunt-malleable

19-195 - 200

Suction curettes, stainless steel, fitting hose connection. (19-201)

Working tips are smoothly rounded off, therefore no danger of perforation.



(suitable for suction biopsy)
19-195: diam 4 mm
19-196: diam 6 mm
19-197: diam 8 mm
19-198: diam 10 mm
19-199: diam 12 mm
19-200: diam 14 mm

19-201

Hose connection inside diam. 11 mm, curved, stainless steel.

(19-195 - 19-200)



OLIVE
 $\varnothing$ 11.15mm

19-201

The new connection between suction curette and hose connection, (19-201) Guarantees quick changing and free turning of the suction curettes on the hose connection.

19-202

Transparent suction hose 1.8 mm long, inside diam 11 mm, outside diam 16.5 mm can be sterilized in an autoclave at 1 bar overpressure (120°C)

19-202

19-203

Strainer for tissues, plastic, for single use only, to be used for histopathological examination.

19-203





19-204 diam 8 mm-19-205 diam 10 mm
BOZEMANN-FRITSCH
 Uterine Flushing Tube
 350 mm



19-207



19-208

- 19-206**
SCHULZE
 Salpingograph, for utero-salpingography
 in wooden case, consisting of:
19-207 Manometer only
19-208 Connection only, without manometer
19-209 Uterine cannula with tape without cone
19-210 Cone only, small
19-211 Cone only, medium
19-212 Cone only, large
19-213 Record syringe, 10CC, luer-lock mount.



19-209



19-210



19-211



19-212

19-213



- 19-214**
SCHULZE
 Salpingograph, equipment like(19-206)
 but without manometer.

