


Obstertrics & Tracheotomy

20



PAGE 246 TO 256

www.humcorp.com



**20-101
THOLE**
for internal and external measurements,
graduated in centimeters



**20-102
DE LEE**
graduated in centimeters
and inches

for internal measurements
between symphysis and
promontory



**20-103
COLLIN**
graduated in
centimeters



**20-104
MARTIN**
graduated in
centimeters



**20-105
SCHNEIDER**
graduated in centimeters

The additional straight arm
can be replaced. The pelvimeter
with straight arm is used for
localisation of the portio vaginalis
uteri and of vaginal metastasis.



**20-106
BREISKY**
graduated in
centimeters

Ovum Curettes, Abortion Scoops



Gyneco

- 20-107 4 mm ϕ
- 20-107A 6 mm ϕ
- 20-107B 8 mm ϕ
- 20-107C 10 mm ϕ
- 20-107D 12 mm ϕ
- 20-107E 14 mm ϕ
- 20-107F 16 mm ϕ




Smellie
20-107G
33 cm / 13"



Wallich
20-107H
42 cm / 16 1/2"



20-107J 12 mm
20-107K 14 mm
20-107L 16 mm
20-107M 20 mm
Blunt



Toothed back



20-108
14.5 mm
Blunt



20-108A
18 mm
Blunt




20-108B
21 mm
Blunt



CUZZI
30 cm / 12"



DANIS
20-108C
24 cm / 9 1/2"






20-114
NAEGELE
355 mm




20-115
NAEGELE
400 mm




KIELLAND
20-115A
42 cm



BARNES-NEVILLE
20-115B
43.5 cm
With traction handle



PIPER
20-115C
44.5 cm



SIMPSON-BRAUN
20-116A
305mm



KIELLAND-LUIKART
20-116B
39.5 cm
Solid blades



MCLEAN
20-116C
39.5 cm / 15 1/2"
Solid blades,
short shanks



BARTON
20-116D
35.5 cm / 14"
One hinged blade

For delivery of the
after coming head in
breach presentation



DEMELIN SUZOR
20-117
32 cm



BOERMA
20-117A
28 cm
For caesarean section




SIMPSON
20-118
23.5 cm
For caesarean section



WRIGLEY
20-118A
25 cm
For caesarean section

For caesarean section

For caesarean section




MURLESS
20-119
25.5 cm
Head extractor
with hinged shaft and sliding collar lock



DEWEY
20-119A
35.5 cm
Obstetrical forceps
with axis traction

Obstetrics Forceps



**20-120A
HERFF
membrane
puncturing
forceps
255 mm**

**20-121A
LOWA
membrane puncturing
forceps
260 mm**

**20-120
GREEN-ARMYTAGE
Caesarean section haemostatic
forceps, for edges of the transverse
incision in the lower segment
220 mm**

**20-121
WILLETT
Scalp Flap forceps
185 mm**

**20-122
GAUSS
Scalp Flap forceps
260 mm**




Fig: 1
20-123

Fig: 2
20-124

Fig: 3
20-125

Fig: 1
20-126

Fig: 2
20-127


Fig: 3
20-128

20-123 - 125
WINTER
290 mm

20-126 - 128
WINTER
290 mm

20-129
SAENGER
275 mm

Obstetrics-Placenta and Ovum Forceps



1/4

1/4

20-132

20-133

253

20-130
SAENGER
285 mm

20-131
SAENGER
265 mm

20-132 - 133
LEWKOWITZ
245 mm



1 jaw fenestrated



**20-135
SMELLIE
Perforator
270mm**




**NAEGELE
20-135A
25 cm**



**20-134
KELLY
320mm**





**20-136
BRAUN
Cranioclast
420mm**



**TARNIER
20-136A
44 cm
Basiotribe**

Tracheal Tubes with inner tubes



20-149 - 158
JACKSON
German Silver

Diam. X length	
20-149 Fig. 0= 4x50 mm	
20-150 Fig. 1= 5x48 mm	
20-151 Fig. 2= 6x58 mm	
20-152 Fig. 3= 7x70 mm	
20-153 Fig. 4= 8x78 mm	
20-154 Fig. 5= 9x84 mm	
20-155 Fig. 6= 10x84 mm	
20-156 Fig. 7= 11x86 mm	
20-157 Fig. 8= 12x88 mm	
20-158 Fig. 9= 13x90 mm	

Diam. X length	
20-159 Fig. 1= 5x48 mm	
20-160 Fig. 2= 6x58 mm	
20-161 Fig. 3= 7x70 mm	
20-162 Fig. 4= 8x78 mm	
20-163 Fig. 5= 9x84 mm	
20-164 Fig. 6= 10x84 mm	
20-165 Fig. 7= 11x86 mm	
20-166 Fig. 8= 12x88 mm	
20-167 Fig. 9= 13x90 mm	
20-168 Fig. 10= 14x90 mm	

20-159 - 168
JACKSON
Stainless Steel

Wound-and Tracheal Hooks and Dilators



20-169
BOSE
165 mm


20-170
BOSE
165 mm

20-171 - 173
180 mm

20-174
FINSEN
50 mm

20-175
FINSEN
70 mm

256



20-176
TROUSSEAU
120 mm

20-177
LABORDE
125 mm

20-178
LABORDE
strong pattern
135 mm

20-179
Cleaning Brush
for tracheal tubes
285 mm