

Plaster Instruments



28

PAGE 393 TO 397

www.humcorp.com

Plaster Knives



with ebony handle

28-101
HOPKINS
200 mm



28-102
ESMARCH
180 mm



REINER
28-103
18.5 cm

Plaster Saw



28-104
ENGEL
160 mm



MOD. NEW
28-105
15 cm



SMITH
(Mod. USA)
28-106
20 cm



KNOWLES
28-107 Straight
28-108 Angled
14 cm



LISTER
28-109 (9cm)
28-110 (14cm)
28-111 (18cm)
28-112 (20cm)



UTILITY
28-113 - 121
180mm

28-113	Black	28-117	Yellow
28-114	Red	28-118	White
28-115	Blue	28-119	Orange
28-116	Green	28-120	Voilet
		28-121	Bright Pink

Plaster Scissors



For small plaster dressings



28-122
160 mm



28-123 - 124
ESMARCH
28-123 200mm
28-124 230mm



28-125
BERGMANN
230 mm



BRUNS
28-126 Plain
23.5 cm
T.C. GOLD
28-126A Serrated
24 cm
Open shanks with pin



BRUNS
28-127
23.5 cm
Closed shanks, serrated



BRUNS
28-128
24 cm
Open shanks with pin, serrated




SEUTIN
28-128A
23 cm

Also To Be Used As
Bandage Scissors


Also To Be Used As
Bandage Scissors



TESSIER
28-129
19 cm
1 serrated cutting edge




TESSIER
28-129A
19.5 cm
1 serrated cutting edge




UNIVERSAL
28-130
20 cm


Bandage Scissors



LORENZ
28-131
23 cm



EXCENTRIC
28-131A 16 cm
28-131B 19 cm



TORNADO
28-131C
24 cm
Undulated cutting edge

Plaster of Paris Shears



28-132A
Spare blade



open jaws

28-133A
Spare blade



open jaws



open jaws


Special features:

The construction of the open jaws avoids tiresome clogging by plaster remnants and allows furthermore an easy cleaning. The shears do not glide off the plaster owing to the deep teething on the jaws.

The knives with hollow grinding have a better cutting ability and therefore less force is required.

The knives are guided by ball bearings.

28-132
STILLE
260 mm



28-133
STILLE
360 mm




STILLE
28-134
24 cm



STILLE-MINI
28-135
Blade with reinforced back
28-136
Blade without reinforced back
20 cm




STILLE



28-137 23 cm
Blade with reinforced back
28-137A 23 cm
Blade without reinforced back
28-137B 26 cm
Blade with reinforced back
28-137C 26 cm
Blade without reinforced back
28-137D 37 cm
Blade with reinforced back
28-137E 37 cm
Blade without reinforced back



WOLFF



28-138
18 cm



28-138A
24 cm



BEESON
28-138B
30 cm




HENNIG
28-139
28 cm



HAGLUND
28-140
28 cm



MOD.USA
28-141
23 cm



SMITH
28-142
23.5 cm

